

Plasmid sequencing policy and procedures for BMD KSQ

Updated 2024 January



Plasmid sequencing policy for BMD KSQ

- **ALL** plasmid sequencing will be done in-house, with only a limited set of exceptions.
- Colleagues are expected to use the in-house procedures outlined in this slide deck.
- Colleagues are reminded to
 - Be responsible: provide the colleagues who are doing the sequencing work with advance notice and all the information they need.
 - Be flexible as in-house procedures are adjusted to workload & technology.
- Library sequencing (e.g. bulk NGS of phage outputs or immune repertoires) is handled separately from the procedures discussed in this slide deck. Contact: Jon McDaniel

Starting in January 2024, all in-house sequencing requests will use new GDBxT system

- Separate requests for full plasmid and amplicon sequencing
- Submit request through GDBxT
 - [Details of GDBxT request system here](#)
 - While the system is in the early roll-out period (Jan-Feb 2024), please check with the appropriate contact person (Hamza Sahil, Ilya Tikh, HTP Nanopore team) to make sure that the correct information is in place.
- Bioinformatics contacts:
 - Request system will notify appropriate bioinformatics scientist.
 - Users with questions about the data they receive should contact Lijian Yu for Nanopore, Tatyana Zamkovaya for PacBio

Procedures for BEP, RNA, TPS groups

- Protein expression constructs (without polyA):
 - Contact: Hamza Sahil
 - Process: Full-plasmid sequencing [Details here](#)
- RNA expression constructs with polyA:
 - For requests of >32 plasmids: PacBio. Contact: Ilya Tikh [Details here](#)
 - For requests of <32 plasmids:
 - Full-plasmid sequencing (Contact: Hamza Sahil [Details here](#))
 - PolyA check (Wyzer): single read, Sanger, raw data, no editing. [Details here](#)
- No need to re-sequence vendor constructs that are already sequence-confirmed

Procedures for **Antibody Discovery/Engineering** groups


- Antibody discovery sequencing (e.g. panels from phage, immune sources)
 - Contact: HTP nanopore team (Shayne Sherry, Joe Bedard, Oliver Ho, Lauren Krause)
 - Process: Amplicon sequencing [Details here](#)
- Protein engineering sequencing (e.g. in-house cloned constructs)
 - Protein expression constructs:
 - Contact: Hamza Sahil
 - Process: Full-plasmid sequencing [Details here](#)
 - RNA expression constructs encoding polyA
 - For requests of >32 plasmids: PacBio. Contact: Ilya Tikh [Details here](#)
 - For requests of <32 plasmids:
 - Full-plasmid sequencing (Contact: Hamza Sahil [Details here](#))
 - PolyA check (Wyzer): single read, Sanger, raw data, no editing. [Details here](#)
- No need to re-sequence vendor constructs that are already sequence-confirmed

Procedures for HTP

- Protein expression constructs (without polyA):
 - Contact: HTP Nanopore team (Shayne Sherry, Joe Bedard, Oliver Ho, Lauren Krause)
 - Process: Amplicon sequencing [Details here](#)
- RNA expression constructs encoding polyA
 - For requests of >32 plasmids: PacBio. Contact: Ilya Tikh [Details here](#)
 - For requests of <32 plasmids:
 - Full-plasmid sequencing (Contact: Hamza Sahil [Details here](#))
 - PolyA check (Wyzer): single read, Sanger, raw data, no editing. [Details here](#)
- No need to re-sequence vendor constructs that are already sequence-confirmed

How to use GDBxT request system – Plasmid Sequencing 1

1. Make sure that there are Plasmid Batches (VBs) for each plasmid you want sequenced
 - a) If not you can create them in Genedata either through the VEC page for your plasmid or by bulk upload
 - b) Be sure to enter a concentration for the VB and set the status to 'Prepped'
2. From a table of VB's (e.g. the [Plasmid Batch Compendium](#)) select the samples of interest and choose 'Create NGS QC Request'



The screenshot displays a web interface for managing plasmid batches. At the top, it shows '93534 Plasmid Batches' with links for 'Create LRs', 'Extended Table', and a settings icon. Below this, a yellow bar indicates '4 rows selected' and provides 'Select all rows' and 'Unselect all' options. A 'Create NGS QC Request' button is also visible. The table below lists four selected rows, each representing a plasmid batch. The columns include ID, Name, Synonym, Created With Production Plate, Verification Status, Prepared By, Preparation Date, Preparation Method, DNA Amount [µg], and DNA Concentration [µg/mL].

<input type="checkbox"/>	ID	Name	Synonym	Created With Production Plate	Verification Status	Prepared By	Preparation Date	Preparation Method	DNA Amount [µg]	DNA Concentration [µg/mL]
<input checked="" type="checkbox"/>	VB-93962	BMD-BEP TRNS-VEC-91200-1		No		Cheryl Lynn Habrukowich	2023-Dec-15	giga	60000.0	2000.0
<input checked="" type="checkbox"/>	VB-93961	BMD-BEP TRNS-VEC-91199-1		No		Cheryl Lynn Habrukowich	2023-Dec-15	giga	60000.0	2000.0
<input checked="" type="checkbox"/>	VB-93960	BMD-BEP TRNS-VEC-91198-1		No		Cheryl Lynn Habrukowich	2023-Dec-15	giga	60000.0	2000.0
<input checked="" type="checkbox"/>	VB-93959	BMD-BEP TRNS-VEC-91197-1		No		Cheryl Lynn Habrukowich	2023-Dec-15	giga	60000.0	2000.0

How to use GDBxT request system – Plasmid Sequencing 2

3. You will be taken to a new page in GDBxT

Create NGS Request

Requestor	Joel Bard
Date Submitted	01/02/2024
Purpose*	Plasmid Verification/QC
Description	<input type="text"/> <small>Maximum 255 characters.</small>
Please provide additional info	<div><input type="text"/> <small>2000 characters left</small></div>
Target Concentration(ug/mL)	20
ReArray Volume(ul)	20
Sequencing Technology Used*	Nanopore
Template Type*	Plasmid
Site Used for NGS*	Please select the option
<input checked="" type="checkbox"/> Please check if Multiplexed	
ELN Link*	<input type="text"/> <small>Maximum 255 characters.</small>
XTP-ID	No matching values found

Nice to have

Nice to have

Default is Nanopore. Change to PacBio if needed.

Required. Choose the site where NGS is being done.

Optional. Paste Signals link.
Could be ELN for design or for
DNA prep. This will be revised...

How to use GDBxT request system – Plasmid Sequencing 3

4. Samples must be submitted in a plate. Further down the request page is a table of the VBs that were included in the request. Click Create new XTPlate and a plate will be created in the database. Be sure to put the assigned sample in each well.
5. Click “Submit NGS Analysis Request”. Emails will be sent to the appropriate people.

XTP-ID

Plate Table

XTP-ID	Number of Wells	Clone Type
		plasmid

1 to 1 of 1 |< < Page 1 of 1 > >|

VB Table

XTP-ID	Well	Plasmid Batch (VB) ID	VEC Alias
		VB-93959	pET28b-NHis_wtRNaseIII_codon optimized
		VB-93960	pET28b-CHis_wtRNaseIII_codon optimized
		VB-93961	pET28b-NHis_E117K_RNaseIII_codon optimized
		VB-93962	pET28b-CHis_E117K_RNaseIII_codon optimized

XTP-ID

Plate Table

XTP-ID	Number of Wells	Clone Type
XTP-48396	4	plasmid

1 to 1 of 1 |< < Page 1 of 1 > >|

VB Table

XTP-ID	Well	Plasmid Batch (VB) ID	VEC Alias
XTP-48396	A1	VB-93959	pET28b-NHis_wtRNaseIII_codon optimized
XTP-48396	B1	VB-93960	pET28b-CHis_wtRNaseIII_codon optimized
XTP-48396	C1	VB-93961	pET28b-NHis_E117K_RNaseIII_codon optimized
XTP-48396	D1	VB-93962	pET28b-CHis_E117K_RNaseIII_codon optimized

Frequency	Percentage
Daily	75%
Not Daily	25%

-
- | Frequency | Percentage |
|-----------|------------|
| Often | 60% |
| Not often | 40% |

Frequency	Percentage
Often	60%
Not often	40%

How to use GDBxT request system – Plasmid Sequencing 5

- Consolidation will create a new NGS request where the person executing the run will create a new XtPlateSet which consolidates all the VBs from the parent NGS requests into a single plate for processing

Submit NGS Analysis RequestCreate New XTPlate

Plate Table

XTP-ID	Number of Wells	Clone Type
XTP-48417	7	plasmid


1 to 1 of 1 < > Page 1 of 1 > >|


VB Table

XTP-ID	Well	Plasmid Batch (VB) ID	Parent NGS ID	VEC Alias
XTP-48417	A1	VB-93160	NGS-4	VEC-83465_NtermGLP1var
XTP-48417	B1	VB-93161	NGS-4	VEC-48637_Nterm_GLP1var
XTP-48417	C1	VB-93162	NGS-4	VEC-83465_GLP1_GCG_Nterm_fusion
XTP-48417	D1	VB-93163	NGS-4	VEC-48637_GLP1_GCG_Nterm_Fusion
XTP-48417	E1	VB-93164	NGS-5	ptt5apopeptide
XTP-48417	F1	VB-93165	NGS-5	ptt5scrambleapopeptide
XTP-48417	G1	VB-93166	NGS-5	pTT5 RAP 536 mlgG2A WT GVHS_C...


How to use GDBxT request system – Amplicon (Screening) Sequencing 1

1. This allows you to request sequencing from Genedata plates.
 - a) Go to the relevant PlateSet within your Screening Campaign.
 - b) Select the PLT IDs that you want included in your Nanopore run.
 - c) Click the Gear Icon and select “Create NGS Request from plates”

 4 Plates

[Get Barcode Names](#) [Plate View](#) [Group](#) 

All rows selected [Unselect all](#)

<input checked="" type="checkbox"/>  ID	Name	Description	Type	Plate Format	Plate Layout
<input checked="" type="checkbox"/> PLT-40033	SC-1430-M-012	p7 of sorted plates	Master Plate	96	All filled
<input checked="" type="checkbox"/> PLT-40032	SC-1430-M-011	p5 of sorted plates	Master Plate	96	All filled
<input checked="" type="checkbox"/> PLT-40031	SC-1430-M-010	p3 of sorted plates	Master Plate	96	All filled
<input checked="" type="checkbox"/> PLT-40030	SC-1430-M-009	p2 of sorted plates	Master Plate	96	All filled

Create GATC Sequence Request File

Create SeqTech Sequence Request File

Plate Sequencing Order Form


Align PLT sequences


Create NGS Request from plates


Print labels

Kenta Shinoda

Caution: the names of the clones must currently begin with SC followed by 4 digits (**SC1234**) in order for software to parse sequences. This may change in the future.

 96 Well Contents **good**



<input type="checkbox"/>  ID	Name	For
<input type="checkbox"/> CL-2871959	SC1426-2H06-2-2	hl

 32 Well Contents **Will not work**



<input type="checkbox"/>  ID	Name
<input type="checkbox"/> CL-2873895	83862(QQ)_W_H33_D_W_W95_QVEC-83910_S_L27E_E_S_L56_E_W_L89_Q-5

How to use GDBxT request system – Amplicon (Screening) Sequencing 2

2. You will be taken to a new page in GDBxT

Create NGS Request

Requestor	Joel Bard
Project	PRJ-551 Alk1 mAb Athero
Screening Campaign	SC-1430 anti-human ALK1 FACS single B cell sorting (ATX-GK mic) - +
Date Submitted	12/21/2023
Purpose*	Please select the option Required
Description	<div>Please select the option Plasmid Verification/QC Plasmid Verification with PolyA Sub-Cloning Insert Verification Screening Other</div>
Please provide additional info	<div></div> <div>2000 characters left</div> <div>Nice to have</div>
Sequencing Technology Used*	Nanopore
Template Type*	Antibody Amplicon
Site Used for NGS*	Please select the option Required
<input checked="" type="checkbox"/> Please check if Multiplexed	
ELN Link*	<div></div> <div>Maximum 255 characters.</div> <div>Required. Any relevant Signals link.</div>

How to use GDBxT request system – Amplicon (Screening) Sequencing 3

3. For each plate set the Clone Type and Primer Set. Choice in the top row autofills down.
- Note that if you're sequencing VH and VL from the same PLT IDs ([hybridoma workflow](#)) you'll have to create two NGS requests which will then be consolidated.

Site Used for NGS* BMD-Cambridge

☒ Please check if Multiplexed

ELN Link* https://prod.snbprd1ykv.perkinelmercloud.com/elements/entity/experiment:f623582e-d7a4-4da0-a59a-f6dc1fa3174c?focus

Maximum 255 characters.

Submit NGS Analysis Request Click Submit and emails will be sent

PLT TABLE

SC PLT-ID	Number of Clones	Clone Type	Primer Set
PLT-40030	96	Please Select...	PhoAphage
PLT-40031	96	Please Select...	Please Select...
PLT-40032	96	VH	pWRIL
PLT-40033	96	VL	PhoAphage
		scFv	PhoAphage
		Peptide	

1 to 4 of 4 | Page 1 of 1

How to use GDBxT request system – Amplicon (Screening) Sequencing 4

4. The full set of NGS requests can be viewed through the [NGS Portal](#)
5. If (**and only if**) you are doing the PCR you need to create a “Consolidated” request even if only running a single request. Check with Sara before doing this as there may be other requests that could be included in the same run.
 - a) Select 1 or more plates that will be included in the same Nanopore run and click Consolidate

Nanopore Requests

Request ID	Consolidated	Purpose	Technology	Site	Status	Requestor	Job Group Name	Project
<input type="checkbox"/> NGS-3	<input type="checkbox"/>	Screening	Nanopore	BMD-Cambridge	New	Joel Bard		551 Alk1 mAb Athero
<input type="checkbox"/> NGS-2	<input type="checkbox"/>	Plasmid Verification/QC	Nanopore	BMD-Cambridge	New	Joel Bard		
<input type="checkbox"/> NGS-1	<input type="checkbox"/>	Plasmid Verification/QC	Nanopore	BMD-Cambridge	New	Joel Bard		

1 to 3 of 3

<

>

Page 1 of 1

Consolidate

Delete

How to use GDBxT request system – Amplicon (Screening) Sequencing 5

- 6) You will see a new NGS Request page.
 - a) The main point of this is to assign the reverse barcodes which distinguish the individual plates in the final results.
 - b) You can assign the plate indices at the bottom of the page.
 - c) Be sure to use the primer plates that correspond to the assigned Primer Set Plate Index for each PLT

Mark as Consolidated

PLT TABLE

SC PLT-ID	Number of Clones	Clone Type	Primer Set	Primer Set Plate Index
PLT-40030	96	scFv	pWRIL	1
PLT-40031	96	scFv	pWRIL	2
PLT-40032	96	scFv	pWRIL	3
PLT-40033	96	scFv	pWRIL	4

1 to 4 of 4

<

>

Page 1 of 1

How to use GDBxT request system – Amplicon (Screening) Sequencing 6

- 7) After Mark as Consolidated has been clicked new buttons will appear: Update Execution And Results
- 8) If you click on Execution and Results you'll be able to download a Deck Layout file indicating which plates should go in which position on the i7 for PCR setup

NGS Execution

Library Prep

Assigned to:

Assign to me

ELN for Library Prep:

Sample Prep Kit:

Update

Download Deck Layout 1

Set Library Prep Complete

Primers for Nanopore Amplicon Sequencing

- For amplicon sequencing, all primers need to be registered in the GDBxT database.
- If you need to use primers that are not available in the database please discuss with Sara Hanscom and Joel Bard **before** starting your PCR.
- The list of available primer sets can be found [here](#). New primer sets must be submitted using the [template](#) (click the Download Template tsv button) from GDBxT.

Procedure for Plasmid Sequencing requests

- Request Submission: Requestors must submit their requests using GDBxT by **10AM on Tuesdays and Thursdays**.
- Plasmid Sample Submission: Requestors must submit purified plasmid samples to **USAX3-610N/004/429 in the DNA drop off fridge**.
- To be included in the NGS run for next-day results, samples need to be **dropped off by 12PM on Tuesdays and Thursdays**.
- Plasmid submission details:
 - **Reference Sequences:** **Confirm the correct reference sequences are uploaded to GeneData before sample drop-off.**
 - **Purification Methods:** All purification methods are acceptable. If using a purification kit with a precipitation step, ensure that any **residual alcohol is removed from the sample**.
 - **Concentration:** Each sample **MUST** have a volume of **100uL and be normalized to 20 ng/uL**.
 - **Tube/Plate Format:** Samples **MUST** be submitted in either **8-strip tubes or a 96-well plate, arranged column-wise starting from position A1**.
 - **Sample Naming:** Each sample **MUST** have the VB-ID as the sample name.

Procedure for Amplicon Sequencing requests

Deadlines may change. Please connect with Shayne, Joe, Oliver, and Lauren as early as possible.

- Contacts: HTP nanopore sequencing team (Shayne Sherry, Joe Bedard, Oliver Ho, Lauren Krause).
 - Individuals from discovery groups may also be trained as superusers.
- For phage library clones:
 - Requestor prepares 100 ul of 1:10 dilution of glycerol stocks in water.
 - Requestor submits request in GDBxT (using PLT numbers) and provides diluted glycerol samples to HTP **by 12PM Tuesdays and Thursdays.**
 - HTP completes plate consolidation in GDBxT, runs PCR and Nanopore prep.
- For pTT5/other vector clones made in HTP:
 - Requestor submits request in GDBxT (using PLT numbers) **by 12PM Tuesdays and Thursdays.**
 - Templates are 1:10 glycerol stock dilutions

Procedure for PacBio requests

- Please contact Ilya at least 1 week prior to expected run date, since PacBio is not run on a regular schedule
- Requestor submits request in GDBxT
- Requestor submits purified plasmid samples to Ilya the day before the sequencing run. Data is usually available in 3-5 business days after sample receipt
- Plasmid submission details:
 - Acceptable purification methods: Miniprep, >20ul of sample
 - Concentration: 30-100ng/ul, standardized across the plate
 - Tube/plate format: 96-well preferred, 8-well strips also ok
 - Sample naming: All samples must have VBs registered in GD

Procedure for Wyzer requests

- ONLY the following request type from Wyzer is permitted under the current budget. This is focused on sequencing just the polyA region of RNA vectors
 - Questions on approved uses of Wyzer sequencing should go to May Tam.
- Order online via Wyzer website
 - Create [Sanger Order](#) → Result type → [Raw Data](#)
- This does NOT include editing or double-stranded reads. Users are expected to address questions on their own by inspection of chromatograms supplied by Wyzer.
 - User submits DNA at 100-200 ng/μL. Most commercial kits should be OK; avoid residual EtOH or salt (some care may be needed w. BioTage prep; consult Hamza, other Nanopore users in TPS,BEP)
 - **Primer design for specific plasmids should read through polyA in single reads. Wyzer has designed primers. Consult Lisa Racie, Nathanael Lintner, Ilya Tikh for details on preferred primers.**
 - [S1 2024](#): users should identify primers that are working well for standard BMD vectors and update this slide deck.
 - VEC ____: Primer name ____ Primer sequence ____ Note: reads from 3' vector sequence through polyA

Note: There is a \$15.00 per order charge (per primer, not per reaction) for design/handling/storing custom primers at Wyzer. Raw reads should be charged \$6 per reaction in addition to this per-order charge. Please review [Wyzer guidelines](#) for sample preparation and submission