


Creating, Submitting, and Viewing Genedata Biosensor Request: A Step-by-Step Guide

Paolo Casas and
Kerry Kelleher

02/09/2024






Overview

- Biosensor KSQ Genedata Integrated Workflow
- Creating, Editing, and Deleting Genedata Biosensor Requests
 - Creating Requests from PPB-IDs (Purified)
 - Creating Requests from Plate Sets (Unpurified: hybridoma, TAP sups, etc.)
 - Modifying Existing Biosensor Requests
- Viewing Genedata Biosensor Assay Results

KSQ Biosensor Genedata Workflow



- Workflow starts and ends with Genedata
- Biosensor request facilitates liquid handling/biosensor instrument set up by providing the ligand/analyte information
- Request system provides an organized and centralized location for data to be stored and accessed

A large, abstract graphic composed of several overlapping, curved, and faceted planes in various shades of blue and purple. The planes are arranged in a way that suggests a complex, three-dimensional structure, possibly representing a molecular model or a data visualization. The lighting creates soft gradients and highlights on the surfaces, giving it a sense of depth and volume.

Creating Biosensor Requests from Purified Material (PPB- IDs)

1. Creating Biosensor Requests from Purification Batch Compendia (PPBs)

Genedata Biologics Home → Compendia → Production → Purification Batches

Genedata Biologics/Bioprocess

Compendia

Reports

Search

Bulk

Reso

Biologics Research Project

580 Projects

ID	Name	Description	Status	Start Date
PRJ-587	Common light chain saturation mutagenesis			2024-Jan-3
PRJ-586	Next Gen pre-humanized semi-synthetic VHH Library			2024-Jan-3
PRJ-585	TREM2/CD98 Bispecific			2024-Jan-3
PRJ-584	Amyloid Beta/CD98 Bispecific	AbetaCD98		2024-Jan-3
PRJ-583	Amyloid Beta	Abeta		2024-Jan-3
PRJ-582	CD98	CD98		2024-Jan-3
PRJ-581	Activin-GDF8-GLP1	ActGDF8		2024-Jan-1
PRJ-580	Targeted IL-12_PharmSci	IL-12	Pha Dev	

Projects

Variable Regions

Inserts

Target Products

Plasmids

Screening

Engineering

Production

Antigens

Antigen Batches

Attachments

Requests

Production Datasets

Production Plate Sets

Production Plates

Plasmid Combinations

Plasmid Batches

Glycerol Stocks

Host Cell Line Batches

Cell Lines

Cell Line Batches

Expression Batches


Purification Batches

Laboratory Results

Actual Products

Aliquot Groups

Aliquot Transactions



Select Filter Panel

92327 Purification Batches

Register Create LRs CoA Extended Table

ID	Name	Synonym	Batch Type	Composition Type
PPB-99471	BMD-BEP STBL_GBT-GIPRGLP1R-0025_PB-11_GG1L		Protein	Single Component
PPB-99470	01302024 huALK5-ECD_fXa_monoFc_His6	CID3641	Protein	Single Component
PPB-99469	01302024 huALK1-ECD_3C_monoFc_His6 (no Kif)	CID3763	Protein	Single Component
PPB-99468	Flu-cHhg3/B_H6 PB1		Protein	Single Component

1. Creating Biosensor Requests from Purification Batch Compendia (PPBs)

Set property to “ID”, operator to “In”, and paste a comma or space separated list as values

Note that if you received an excel file with the PPB-IDs listed in the column you can simply select the cells in the column and copy/paste directly into the value field.

Filter columns

Choose filter...

Discard

Clear

All - of the following are true

	Operator	Value(s)		
ID	In	PPB-92440 PPB-92441 PPB-92442 PP	x ?	+ -
Choose property	Operator	Value(s)		+ -

Show filter expression

Save as New

Delete

Close

Apply

1. Creating Biosensor Requests from Purification Batch Compendia (PPBs)

Pressing “apply” filters the purification batches, check all PPBs


9 (out of 92327) Purification Batches Register Create LRs CoA Extended Table ⚙ ...

All rows selected ✕

<input checked="" type="checkbox"/> ▾	ID	Name	Batch Type	Composition Type	Use as Antigen Material	Concentration
<input checked="" type="checkbox"/>	PPB-92448	BMD-CD27-41BB-4004_Glu_4W	Protein	Complex	No	5.07
<input checked="" type="checkbox"/>	PPB-92447	BMD-CD27-41BB-4004_Glu_2W	Protein	Complex	No	5.02

Indicates filter is applied

Gear Icon will allow creation of Biosensor Request from selected PPBs



ister Create LRs CoA Extended Table ⚙

Register PF Number for Purification Batch

Add PPBs to Existing Analysis Request Set

Alignment and properties export

Request Biosensor Analysis From PPBs

Request Formulation for RNA

Request Mass Spec Analysis

Request Molecular Assessment

Request mRNA Molecular Assessment

Export PPBs to Excel

1. Creating Biosensor Requests from Purification Batch Compendia (PPBs)

Enter Analysis Request Set Details within GDBxT

Analysis Request Set Details

Project * PRJ-530 | GIPR

Date Submitted * 02/09/2024 Enter date of analysis request set submission

Requesting Scientist * casasj07 Enter your Genedata username

BMD Team Leader James Apgar

Biology Team Leader Please enter a name

Purpose * purpose goes [here](#) Enter the reason for requesting this analysis
983 characters left

Analysis * Off-rate

Analysis Group * MDPP - Biosensor

Project Lead From Group * Kerry Kelleher

Comments comments can go [here](#) Enter any additional instructions for the analysis group
235 characters left

- * Fields are required with many being prepopulated
- Enter the **Purpose** of the request. Although not required, please provide the binding partner species, concentration, TP Name, PPB ID,
- Select the the **Analysis** type, (K_D Values, Off-rate,...)
- Select the **Project Lead from Group**.
- Include **Comments** that will be helpful for the Biosensor scientist to complete the request
- Information regarding which instrument was used is not required at this time, but will later be included as part of the laboratory results

Note: A request is a pairing of two proteins (example 1 Ab and 1 Ag), whereas a Request Set may include multiple pairings (example 5 Abs against 1 Ag).

1. Creating Biosensor Requests from Purification Batch Compendia (PPBs)

Complete the Request by Selecting the Protein Purification Batch (PPB)

Submit the Request by Selecting the 'Create Request Set' Button

The screenshot displays a web interface for creating biosensor requests. It features two rows of protein purification batches (PPBs). Each row includes a checkbox, a protein ID (e.g., TPP-110792, GBT-GIPR-3578), and a dropdown menu showing the selected PPB (e.g., PPB-91681 | GBT-HTP_GBT-GIPR-3578_PB-1). Below the rows is a green button with a '+' icon and a red arrow pointing to a blue 'Create Request Set' button. A red arrow also points to the dropdown menu in the first row.

Checkbox	Protein ID	PPB Selection	Buttons
<input type="checkbox"/>	13 TPP-110792 GBT-GIPR-3578	0.92 mg/mL PPB-91681 GBT-HTP_GBT-GIPR-3578_PB-1	Add PPB Reduce to one PPB
<input type="checkbox"/>	14 TPP-110793 GBT-GIPR-3579	1.38 mg/mL PPB-91682 GBT-HTP_GBT-GIPR-3579_PB-1	Add PPB Reduce to one PPB

Buttons: +, Create Request Set, Remove Proteins

By requesting through PPBs in the purification batch compendia, the exact TPP-PPB hierarchy is automatically supplied for the user

Example of GDBxT Message Displayed after Selecting the 'Create Request Set' Button

Analysis Request Set 2089 Successfully Submitted To Biosensor Group

1. Creating Biosensor Requests from Purification Batch Compendia (PPBs)

Occasionally, the user will still need to select the specific PPB from the TPP using this method (ex. child PPBs with the same parental PPB ID/TPP for forced degradation analysis or bispecifics)

<input type="checkbox"/>	TPP ID	TPP Name	Protein Purification Batch
<input type="checkbox"/>	1	TPP-105941	BMD-CD27-41BB-4004

+

Create Request Set

Remove Proteins

Cancel

Add PPB

Reduce to one PPB

1.0 mg/mL | PPB-86794 | BMD-LJ_BMD-CD27-41BB-4004-1

4.61 mg/mL | PPB-92441 | BMD-CD27-41BB-4004_Trīs_2W

121.7 mg/mL | PPB-91981 | BMD-BEP STBL_BMD-CD27-41BB-4004-4

121.7 mg/mL | PPB-91981 | BMD-BEP STBL_BMD-CD27-41BB-4004-4

4.8 mg/mL | PPB-89873 | BMD-BEP STBL_BMD-CD27-41BB-4004-1

4.88 mg/mL | PPB-92442 | BMD-CD27-41BB-4004_Trīs_4W

4.91 mg/mL | PPB-92440 | BMD-CD27-41BB-4004_Trīs_T0

4.92 mg/mL | PPB-92443 | BMD-CD27-41BB-4004_His_T0


4.95 mg/mL | PPB-92444 | BMD-CD27-41BB-4004_His_2W

5.0 mg/mL | PPB-92446 | BMD-CD27-41BB-4004_Glu_T0

5.02 mg/mL | PPB-92447 | BMD-CD27-41BB-4004_Glu_2W


(Add or remove a PPB with the same TPP ID)

(Select the PPB for Biosensor analysis. The PPB ID and/or concentration (mg/mL) will help to identify the specific Protein Purification Batch for Biosensor analysis)



2. Creating Biosensor Requests from HTP Protein Production Request Sets

Select HTP Request Sets within GDBxT




Select the HTP Request Set Name

GBT-HTP Home Page

Request Sets							
Active Requests							
Expression Tasks							
Purification Tasks							
Available CM							
Request Set Platesets							
Export Basic							
<input type="checkbox"/>	RS ID	Date Submitted	Request Set Type	Project	AT Code	Formats	Rationale
<input type="checkbox"/>	RS-23428	12/06/2021	-	TGFbeta	AR-700001320	One Armed IgG	RS of 47 + 1 Control TGFB OAA Microscale: data acquisition for Monica
<input checked="" type="checkbox"/>	RS-23427	12/06/2021	-	TGFbeta	AR-700001320	One Armed IgG	RS of 321 TGFB OAA Microscale: data acquisition for Monica
<input type="checkbox"/>	RS-23414	11/29/2021	-	CD3	HD-EX00033	One Armed IgG	screen for super low affinity aCD3 variants

2. Creating Biosensor Requests from HTP Protein Production Request Sets

Select All TPPs under Protein Purification Batches Section (usually found near the bottom of the RS page)



<input checked="" type="checkbox"/>	TPP-ID	TPP Name	PEB-ID	PEB Name	PPB-ID	PPB-Name
<input checked="" type="checkbox"/>	TPP-77980	GBT-TGFB-0740	PEB-49581	GBT-HTP_GBT-TGFB-0740_EB-1	PPB-58756	GBT-HTP_GBT-TGFB-0740_PB-1
<input checked="" type="checkbox"/>	TPP-77979	GBT-TGFB-0739	PEB-49582	GBT-HTP_GBT-TGFB-0739_EB-1	PPB-58757	GBT-HTP_GBT-TGFB-0739_PB-1
<input checked="" type="checkbox"/>	TPP-77978	GBT-TGFB-0738	PEB-49583	GBT-HTP_GBT-TGFB-0738_EB-1	PPB-58758	GBT-HTP_GBT-TGFB-0738_PB-1

Select the 'Request Biosensor Analysis' Button




[Add SEC Results](#) [Add Lab Results](#) [Download File for Tube Labeller](#) [Request Molecular Assessment](#) [Request Biosensor Analysis](#)


Most HTP request sets tend to have a 1:1 relationship between TPP and PPB-ID you will not need to manually select the PPB


3. Creating Biosensor Requests from Target Products (TP)


Genedata Biologics Home Page Search for the Project

 Genedata Biologics/Bioprocess


Compendia ▾ ...


FN14 


 ▾

 ▾

Biologics Research Project Compendium





 450 Projects


Create 






Select the Genedata Project Name



Search Results for 'FN14'

 Open all sections

 Close all sections



 1 Project



    






	ID	Name	Synonym	Target	Description	Status	Technical Project Lead	Internal Project ID	User Group ID	Group	Registered By	Registered Via	Modified By	Modified Via
	PRJ-343	Fn14 Monovalent Antibody Cachexia		FN14		LD	James Apgar	HD-OB01219		Research	May Tam		May Tam	

3. Creating Biosensor Requests from Target Products (TP)

Select the Target Products Filter Icon

 1044 Target Products 

Plasmid Report Create Request Protein Analyzer Create  

3 row(s) selected Select all rows Unselect all										   	
<input type="checkbox"/>  ID	Name	Sequence Locked	Type	Format	Chain Multiplicity	Isotype	Antibody Library ID	Variable Region ID	Constant Reg		
<input type="checkbox"/> TPP-77336	GBT-FN14-2170	No	Antibody	scFv	1	Kappa	LIB-1	VR-13805			
<input type="checkbox"/> TPP-75594	GBT-FN14-2169	No	Antibody	IgG	1, 1	Kappa, IgG1	LIB-1		CON-6, CON		
<input type="checkbox"/> TPP-75593	GBT-FN14-2168	No	Antibody	IgG	1, 1	Kappa, IgG1	LIB-1		CON-6, CON		
<input type="checkbox"/> TPP-75592	GBT-FN14-2167	No	Antibody	IgG	1, 1	Kappa, IgG1	LIB-1		CON-6, CON		
<input type="checkbox"/> TPP-75591	GBT-FN14-2166	No	Antibody	IgG	1, 1	Kappa, IgG1	LIB-1		CON-6, CON		
<input type="checkbox"/> TPP-75533	GBT-FN14-2165	No	Antibody	IgG	1, 1	Kappa, IgG1	LIB-1		CON-6, CON		

Frequency	Percentage
Often	50%
Not often	50%

Filter for the TD (biotherapeutic) Name

Filter columns

Choose filter...

Discard

Clear

Any - of the following are true

(All, Any)

	Operator	Value(s)		
Name	Contains	347	x 🔍	+ - (Add, remove)
Name	Contains	391	x 🔍	+ -
Name	Contains	402	x 🔍	+ -
Name	Contains	0001	x 🔍	+ -

Show filter expression

Save as new

Delete

Close

Apply


Please note: The aim of this search is to identify the biotherapeutics (example Ab), for Biosensor analysis. It is not necessary to select the other binding partner (example Ag) when creating a Genedata Biosensor request.




3. Creating Biosensor Requests from Target Products (TP)

Select the Target Product Proteins for Biosensor Analysis

Using the Gear Icon, Select 'Request Biosensor Analysis' from the Dropdown List

 4 (out of 1044) Target Products

Plasmid ReportCreate RequestProtein AnalyzerCreate 

3 row(s) selectedSelect all rowsUnselect all

<input type="checkbox"/>	ID	Name	Synonym	Description	Type
<input type="checkbox"/>	TPP-61048	GBT-FN14-0402	PF-07860725	H clone: GBT-FN14-0092_heavyGraft2 L clone: GBT-FN14-0092_lightGraft5 MnM FN14 Hum. VH-VLIgG	Antibody
<input checked="" type="checkbox"/>	TPP-61037	GBT-FN14-0391	PF-07860724	H clone: GBT-FN14-0084_heavyGraft4 L clone: GBT-FN14-0084_lightGraft5 MnM FN14 Hum. VH-VLIgG	Antibody
<input checked="" type="checkbox"/>	TPP-57989	GBT-FN14-0347	PF-07860726	H clone: sc809_WBP3840-P9R41-1E6 L clone: sc809_WBP3840-P9R41-1E6 MnM Phage Campaign Wuxi FN14IgG	Antibody
<input checked="" type="checkbox"/>	TPP-43337	GBT-FN14-0001	PF-06756852	H clone: huP4A8 L clone: huP4A8 MnM P4A8IgG	Antibody

Show empty columns

Rows per page: 201-4 of 4<<<>>>

Create Cell Line Bulk Excel File

Create Gene Synthesis Order from Vectors

Request Protein Production

Add Samples to Existing Analysis Request Set


Align sequences

Alignment and properties export

Request Biosensor Analysis

Request Biosensor Analysis




Since only the TPP was selected, the user will have to manually select the specific PPBs for analysis

An abstract 3D graphic composed of several blue, curved, and faceted planes that overlap and curve across the top and right side of the slide, creating a sense of depth and movement.

Creating Biosensor Requests from Unpurified Material (Hybridoma, TAP sups, etc.)

4. Creating Biosensor Requests from Plate Sets (unpurified material)

Search for plate set ID in Genedata Biologics, you can also navigate to the plate set from the screening campaign page



SC-1250 > Plate Set PS-5863 (Confirm plates set (284))

Copy ▾ Re-Array ▾ Reformat ▾ Re-Streak Co-Transfect ▾ Worklist Control File ▾ Assign Antigen Create ACBs ⚙️ ...

Details

ID	PS-5863
Name	Confirm plates set (284)
Screening Stage	Confirmatory
Parent Plate Set	PS-5862 (1st fusion group (156 plates))
Parent Clone Pool	CP-6191 (TSR1450 Hybridoma 120 plates), CP-6192 (TSR1450 Hybridoma 36 plates)
Plate Type	Compression Plate
Plate Format	96
Plate Layout	GBT 96 Well 8 Controls

Show empty properties

Last Modified Date User Groups (no restrictions)

3 Child Plate Sets

Once on the Plate Set page, select “Request Biosensor Analysis for PlateSet” under the top gear icon

Notes:

- 1) This is the only gear icon on the page which has this option
- 2) Selecting the specific clones under the well contents section will not affect which clones are brought over to the request set (you can specifically select clones in the next step in GDBxT)

4. Creating Biosensor Requests from Plate Sets (unpurified material)

Enter ARS Details, note that a new screening campaign field is added

Analysis Request Set Details

Project *

PRJ-530 | GIPR

Screening Campaigns *

SC-1250 | GIPR(ECD)-mRNA immunized Alloy

Analysis Type *

Kd Values

Date Submitted *

01/31/2024

Enter date of analysis request set submission

Requesting By *

casasj07

Enter your Genedata username

BMD Team Leader

James Apgar

Biology Team Leader

James Apgar

Purpose *

Enter the reason for requesting this analysis

1000 characters left

Analysis Group *

MDPP - Biosensor


Assigned to *

Kerry Kelleher

Comments

Enter any additional instructions for the analysis group

255 characters left



Analysis Requests

Search

<input type="checkbox"/>	Plate	Well Address	Clone Id	Clone Name
<input checked="" type="checkbox"/>	PLT-36108	A05	CL-2586949	SC1250-12A09
<input checked="" type="checkbox"/>	PLT-36108	D05	CL-2587029	SC1250-13A01
<input checked="" type="checkbox"/>	PLT-36108	E07	CL-2587608	SC1250-19E08
<input type="checkbox"/>	PLT-36108	A01	CL-2586014	SC1250-1D09
<input type="checkbox"/>	PLT-36108	B08	CL-2587737	SC1250-21A05
<input checked="" type="checkbox"/>	PLT-36108	D04	CL-2586761	SC1250-9H08
<input checked="" type="checkbox"/>	PLT-36108	A09	CL-2587922	SC1250-23B03

Scrolling below the ARS details, you can select the clones you want added to the ARS and press create ARS at the bottom of the page. Note that for plate sets, the unique identifiers in the database are the plate ID and well address.

An abstract 3D graphic composed of several blue, curved, and faceted planes that overlap and curve across the top and right side of the slide, creating a sense of depth and movement.

Modify Existing Biosensor Requests

Editing Biosensor Request Set Details

Select the 'Edit Details' Button on the Analysis Request Set Details Page


After Editing, Select the Check Mark Button, then the 'Save Changes' Button

Analysis Request Set Details

ID	ARS-3117
Project Name	GIPR
Project Activity Tracker Code	HP-C567A000
Analysis Group	MDPP - Biosensor
Date Submitted	10-10-2023
Date Completed	
Request Set Status	Complete
Purpose	Testing out VBA automation for Carterra LSA based results. Using 160 PPB-IDs sent to LJ for Biacore vs. Carterra comparison.
Number of Samples	161
Number of Requests	161
Requesting Scientist	Jorge Paolo Casas
Assigned Scientist	Jorge Paolo Casas
BMD Team Leader	James Apgar
Biology Team Leader	
Production Dataset:	https://gdbio.pfizer.com:8081/Biologics/show/PDS-7898
Comments	

Edit Details

Edit Request Set



ID	ARS-3117
Project Name	GIPR
Project Activity Tracker Code	HP-C567A000
Analysis Group	MDPP - Biosensor
Date Submitted	10-10-2023
Date Completed	Empty
Request Set Status	Complete
Purpose	Testing out VBA automation for Carterra LSA based results. Using 160 PPB-IDs sent to LJ for Biacore vs. Carterra comparison.
Number of Samples	161
Number of Requests	161
Requesting Scientist	Jorge Paolo Casas
Assigned Scientist	Jorge Paolo Casas
BMD Team Leader	James Apgar
Biology Team Leader	Empty
Production Dataset:	https://gdbio.pfizer.com:8081/Biologics/show/PDS-7898
Comments	Write a new comment

(Adding more details to the **Comments** field)

This button will allow for Biosensor group addition of binding partner (s) and control (s), covered in next slide



✓

✕

Edit Request Set

Save Changes

Discard Changes

Adding Binding Partner/Antigen

Select the 'Edit Request Set' Button on the Analysis Request Set Details Page

Select checkbox corresponding to correct antigen, PPB-ID and concentration

have been provided to facilitate selection

Analysis Request Set Details

ID

ARS-3117

Project Name

GIPR

Project Activity Tracker Code

HP-C567A000

Analysis Group

MDPP - Biosensor

Date Submitted

10-10-2023

Date Completed

Request Set Status

Complete

Purpose

Testing out VBA automation for Carterra LSA based results. Using 160 PPB-IDs sent to LJ for Biacore vs. C arterra comparison.

Number of Samples

161

Number of Requests

161

Requesting Scientist

Jorge Paolo Casas

Assigned Scientist

Jorge Paolo Casas

BMD Team Leader

James Apgar

Biology Team Leader

Production Dataset:

<https://gdbio.pfizer.com:8081/Biologics/show/PDS-7898>

Comments

Edit Details

Edit Request Set

Antigen Batch *

	Antigen	Description	Antigen Batch	Batch Concentration	Contact
<input checked="" type="checkbox"/>	cyGIPR_ECD_3c_Avi_TEV_CH23LSFc_His6	cyno GIPR ECD	PPB-81377 cyGIPR_ECD_3c	3.04 mg/mL	Chris Connors
<input type="checkbox"/>	GIPR2761VHVL	IgVHss_GIPR2761VHVL_TEV_monoFc_His6	PPB-91340 GIPR2761VHVL_TEV_monoFc_His6 PB1	8.8 mg/mL	David Hokanson
<input type="checkbox"/>	GIPR1362ECDFab	IgVHss_GIPR-ECD_GIPR1362FabHC_His6 + IgVHss_GIPR1362LC	PPB-92611 GIPR1362ECDFab	18.2 mg/mL	Emily Longo
<input type="checkbox"/>	GIPR1362ECDFab	IgVHss_GIPR-ECD_GIPR1362FabHC_His6 + IgVHss_GIPR1362LC	PPB-92747 GIPRecd_GIPR1362fab fusion PB2	4.3 mg/mL	David Hokanson
<input type="checkbox"/>	GIP Peptide	Peptide Custom Synthesized	PPB-80466 GIP Peptide	unavailable	Susan Benard
<input type="checkbox"/>	GIPR0107ECDFab	IgVHss_GIPR-ECD_GIPR0107FabHC_His6+IgVHss_GIPR0107LC	PPB-92609 GIPR0107ECDFab	17.4 mg/mL	Emily Longo

+

Update Request Set

Remove Proteins

Cancel

Scroll to bottom of page and press “update request set”, now the requests on the main ARS page should display both the ligand and the analyte to complete the pairing

Troubleshooting Adding Binding Partners or Controls

Antigen Batch *	Antigen	Description	Antigen Batch	Batch Concentration	Contact
<input checked="" type="checkbox"/>	cyGIPR_ECD_3c_Avi_TEV_CH23LSFc_His6	cyno GIPR ECD	PPB-81377 cyGIPR_ECD_3c	3.04 mg/mL	Chris Connors
<input type="checkbox"/>	GIPR2761VHVL	IgVHss_GIPR2761VHVL_TEV_monoFc_His6	PPB-91340 GIPR2761VHVL_TEV_monoFc_His6 PB1	8.8 mg/mL	David Hokanson
<input type="checkbox"/>	GIPR1362ECDFab	IgVHss_GIPR-ECD_GIPR1362FabHC_His6 + IgVHss_GIPR1362LC	PPB-92611 GIPR1362ECDFab	18.2 mg/mL	Emily Longo
<input type="checkbox"/>	GIPR1362ECDFab	IgVHss_GIPR-ECD_GIPR1362FabHC_His6 + IgVHss_GIPR1362LC	PPB-92747 GIPRecd_GIPR1362fab fusion PB2	4.3 mg/mL	David Hokanson
<input type="checkbox"/>	GIP Peptide	Peptide Custom Synthesized	PPB-80466 GIP Peptide	unavailable	Susan Benard
<input type="checkbox"/>	GIPR0107ECDFab	IgVHss_GIPR-ECD_GIPR0107FabHC_His6+IgVHss_GIPR0107LC	PPB-92609 GIPR0107ECDFab	17.4 mg/mL	Emily Longo

If you do not see the desired antigen in the table, ensure that the “Use as Antigen Material” field is marked as “Yes” in the PPB page on Genedata for the antigen to populate in table

If you need to add a control sample from a different project to the ARS, you will need to link it’s TPP to your current project

PRJ-517 > Target Product TPP-115651 (RBD XBB.1.5 Biotin) Edit

Add Sequence ▾

Create Related Antigen

Link to Project

Link to Antibody Clone

Unli

Details

ID

TPP-115651

Name

RBD XBB.1.5 Biotin

Details

ID

PPB-93410

Name

RBD XBB.1.5 Biotin

Synonym

Biotinylated RBD XBB.1.5 Omicron

Description

SARS-CoV-2 Spike RBD His Avitag Biotin

Use as Antigen Material

Yes

Show empty properties

Deleting Requests within a Request Set Details

Select the Request (s) to Delete

Select the 'Remove Proteins' Button then the 'Create Request Set' Button

The screenshot shows the 'Request Set Details' interface. A red arrow points to the checkbox for request 161 (TPP-104517, GBT-GIPR-2629). The interface includes buttons for 'Add PPB' and 'Reduce to one PPB'. A dropdown menu is open, showing two options: '0.69 mg/mL | PPB-86301 | GBT-HTf' and '0.69 mg/mL | PPB-86301 | GBT-HTP_GBT-GIPR-2629_PB-1'. Below the dropdown, there are checkboxes for various analysis types: 'Off-rate', 'High Low Affinity', 'Kd Values' (checked), 'Competition', 'CFCA', 'Other', 'Percentage Activity', 'Steady State Affinity', 'Single Cycle Kinetics', 'Octet Quantification', and 'Kinexa'. A blue arrow points from the 'Kd Values' checkbox to a second screenshot. The second screenshot shows a green '+' button, a red arrow pointing to the 'Update Request Set' button, and a red 'Remove Proteins' button. Below these buttons are 'Cancel' and 'Update Request Set' buttons.

Note that users have the option to select multiple analysis types using the checkboxes on the right for a particular TPP. However, we have found that keeping limiting an ARS to a single assay type to a single type works best. If you would like to run additional assays on the same samples, then you can create new ARS.

Deleting a Request Set within Request Set Details

To Delete a Submitted Request, Select the ‘Delete Analysis Request Set’ Button

Analysis Requests

☐ Reverse RequestSet

Search

<input checked="" type="checkbox"/>	AR ID	TPP ID	Ligand TPP Name	Ligand TPP Synonym	TPP Format	Ana
<input checked="" type="checkbox"/>	AR-13257	TPP-43337	GBT-FN14-0001	PF-06756852	IgG	Kd
<input checked="" type="checkbox"/>	AR-13258	TPP-57989	GBT-FN14-0347	PF-07860726	IgG	Kd
<input checked="" type="checkbox"/>	AR-13259	TPP-61037	GBT-FN14-0391	PF-07860724	IgG	Kd

Edit requests

Delete Selected Analysis Requests

Export Selected Requests For Liquid Handling


Upload Excel To Lab Results

Upload PowerPoint

Delete lab results from Genedata

The “Reverse RequestSet” button is provided in the event you plan on reversing the assay (i.e. on the Biacore you would capture the analyte and flow the ligand over the surface). Pressing this button will cause the biotherapeutic samples to be listed under the analyte column and the antigen/binding partner under the ligand column. However, we have experienced complications uploading our data when this button is pressed so we tend to leave this button unchecked and reverse the samples in our custom macro.


Delete Analysis Request Set

A large, abstract graphic on the right side of the slide. It consists of several overlapping, curved, and faceted shapes in various shades of blue and purple, creating a sense of depth and movement. The shapes appear to be part of a larger, complex structure, possibly representing a molecular model or a data visualization.

Viewing Genedata Biosensor Assay Results

Viewing Biacore Kinetic Assay Results on the Request Set Details Page

Select the appropriate biosensor requests page from the biosensor analysis request sets drop down menu in GDBxT



Search for the ARS with Text (exact match) and Date Field Parameters

Select the ARS ID

All Analysis Request Sets For MDPP - Biosensor

Analysis Request Set Index

Abs 347

Between Dates:
Start: 08/01/2021End: 08/31/2021




<input type="checkbox"/>	ID	Group	Project	# Requests	Submitted	Requested By	Assigned Scientist	Purpose
<input type="checkbox"/>	ARS-1917	MDPP - Biosensor	Fn14 Monovalent Antibody Cachexia	5	08/2021/26	Kerry Kelleher	Kerry Kelleher	Epitope bin with sandwich assay and human cleaved FN14 (11) lead Abs 347, 391, 402, and positive control Abs 0001 and 0007.
<input type="checkbox"/>	ARS-1909	MDPP - Biosensor	Fn14 Monovalent Antibody Cachexia	5	08/2021/23	Kerry Kelleher	Kerry Kelleher	Determine KD values for cyno FN14 vs. Abs 347, 391, 402 and positive control Abs 0001 and 0007. Ab0007 has mouse Fc.
<input checked="" type="checkbox"/>	ARS-1876	MDPP - Biosensor	Fn14 Monovalent Antibody Cachexia	10	08/2021/03	Kerry Kelleher	-	Determine KD values of human and mouse FN14 vs Abs 347, 391, 402 and pos controls 0001 and 0007.

Viewing Biacore Kinetic Assay Results on the Request Set Details Page

Analysis Request Set 1876

Analysis Request Set Details

Results Biacore KD

Id	Name	Analyzed Entity ID	Measurement Type	Grouping Criterion	--Timestamp--	Orientation	ka (1/Ms)	kd (1/s)	KD (pM)	t1/2 (s)	Rmax (RU)	Chi² (RU²)	Chi²/Rmax (%)	Trace	Binding Partner
LR-172617	BMD-BEP TRNS_GBT-FN14-0347_PB-3 BMD - TPS_GBT-FN14-0011_PB-1_KD-A-1	PPB-52093	Biacore KD	BMD - TPS_GBT-FN14-0011_PB-1	2021-08-13 00:00:00.0	TPP as ligand	5.18E5	2.26E-3	4360.00	3.07E2	31.90	4.61E-1	1.45		PPB-32775
LR-172618	BMD-BEP TRNS_GBT-FN14-0347_PB-3 BMD - TPS_GBT-FN14-0011_PB-1_KD-C-1	PPB-52093	Biacore KD	BMD - TPS_GBT-FN14-0011_PB-1	2021-08-13 00:00:00.0	TPP as ligand	5.54E5	3.36E-3	6060.00	2.06E2	33.50	3.99E-2	.12		PPB-32775
LR-172619	BMD-BEP TRNS_GBT-FN14-0391_PB-3 BMD - TPS_GBT-FN14-0011_PB-1_KD-A-1	PPB-52094	Biacore KD	BMD - TPS_GBT-FN14-0011_PB-1	2021-08-13 00:00:00.0	TPP as ligand	3.58E7	9.79E-3	274.00	7.08E1	14.60	1.47E-1	1.00		PPB-32775

Viewing Biacore Kinetic Assay Results from the Production Dataset (PDS) Page

Select the Production Dataset link on the Biosensor Request Set Details Page

Analysis Request Set Details

ID

ARS-1876

Project Name

Fn14 Monovalent Antibody Cachexia

Project Activity Tracker Code

HD-OB01219

Analysis Group

MDPP - Biosensor

Date Submitted

08-03-2021

Date Completed

Request Set Status

Complete

Purpose

Determine KD values of human and mouse FN14 vs Abs 347, 391, 402 and pos controls 0001 and 0007.

Number of Samples

5

Number of Requests

10

Requesting Scientist

Kerry Kelleher

Assigned Scientist

Kerry Kelleher

BMD Team Leader


James Appar

Biology Team Leader

Production Dataset:

<http://gdbio.pfizer.com:8080/Biology/show/PDS-4283>

Select the Extended Table Button within the PPB Table

 3 Purification Batches

Create/Link APPs

Create LRs

Create PEBs

Pool


Complex





Link


Unlink


Edit

Extended Table









<input type="checkbox"/> 	ID	Name	Description	Batch Type	Concentration [mg/mL]	Group
<input type="checkbox"/>	PPB-52095	BMD-BEP TRNS_GBT-FN14-0402_PB-3	Lead Clones for analysis	Protein	4.4	Research
<input type="checkbox"/>	PPB-52094	BMD-BEP TRNS_GBT-FN14-0391_PB-3	Lead Clones for analysis	Protein	4.23	Research
<input type="checkbox"/>	PPB-52093	BMD-BEP TRNS_GBT-FN14-0347_PB-3	Lead Clones for analysis	Protein	2.92	Research

 Show empty columns

Rows per page: 20

1-3 of 3







Viewing Biacore Kinetic Assay Results from the Production Dataset (PDS) Page

Select the Sensorgram for Zoom View

Laboratory Results for 3 Purification Batches

(Aggregate Lab Results for PPB in addition to Biacore Assay)

Purification Batch			Analytical HIC : 1	Biacore KD : BMD - TPS_GBT-FN14-0011_PB-1				
<input type="checkbox"/> ID	Name	Target Product I...	Percent Main Peak [%]	Orientation (list)	ka (1/Ms) [1/(...	kd (1/s) [...	KD (nM) [n...	Trace (list)
<input type="checkbox"/> PPB-52095	BMD-BEP TRNS_GBT- FN14- 0402_PB-3	TPP-61048	100.0	TPP as ligand, TPP as ligand	1.54E7	5.2E-3	3.42E-1	
<input type="checkbox"/> PPB-52094	BMD-BEP TRNS_GBT- FN14- 0391_PB-3	TPP-61037	100.0	TPP as ligand, TPP as ligand	4.48E7	1.12E-2	2.56E-1	



Viewing Biacore PowerPoint Report Results from the Production Dataset (PDS) Page

Select the 'Attachment' Button



Select the Reference to Document File Name to View the Biosensor PowerPoint Report



<input type="checkbox"/> ID	Reference to Document	Registered By	Registration Date
<input type="checkbox"/> ATT-27561	Biosensor/Attachments/20210813-17TitrationhumuCleavedFN14vsAbs347391402PosCon00010007CM58K37C.pptx	Kerry Kelleher	2021-Aug-22 13:03:45 -0400

Biosensor PowerPoint Report

Biacore Characterization of Human and Mouse FN14 and Abs 347, 391, 402, and Positive Control Abs 0001 and 0007

Kerry Kelleher, Justin Cohen, and Laura Lin

08/22/2021

Pfizer BioMedicine Design

Pfizer Confidential - Internal use only