RE: GeneData Error for new VSIG4 sequences Subject:

Date: Thursday, November 16, 2023 at 2:39:27 PM Eastern Standard Time

From: Sausen, Christopher

To: Bard, Joel

Jackobek, Ryan Andrew, Francis, Christopher Alan, Bitzas, Gabrielle, Kerns, Kelvin

Attachments: image002.png, image005.png, image006.png, image007.png

Hi Joel,

Thanks so much, that worked perfectly. Easy fix!

Best, Chris

## **Christopher Sausen**

Senior Scientist

#### **BioMedicine Design**

**Biotherapeutic Production** 

O: 1-978-247-1374 M: 1-218-851-7391

Christopher.sausen@pfizer.com



From: Bard, Joel < <a href="mailto:Joel.Bard@pfizer.com">Joel.Bard@pfizer.com</a>> Sent: Thursday, November 16, 2023 1:31 PM

**To:** Sausen, Christopher < <a href="mailto:Christopher.Sausen@pfizer.com">Christopher.Sausen@pfizer.com</a>>

Cc: Jackobek, Ryan Andrew < Ryan.Jackobek@pfizer.com >; Francis, Christopher Alan < Christopher.Francis@pfizer.com >; Bitzas, Gabrielle < Gabrielle.Bitzas@pfizer.com >;


Kerns, Kelvin < Kelvin. Kerns@pfizer.com>

Subject: RE: GeneData Error for new VSIG4 sequences

Hi Chris and Ryan-

The way to change complexation is to use the Replace button on the sequence tab on the TPP page and choose the Copy Paste option. That will let you edit the complexation. You won't need to do anything to the sequence, just change complexation and Finish.





From: Sausen, Christopher < <a href="mailto:Christopher.Sausen@pfizer.com">Christopher.Sausen@pfizer.com</a>>

Sent: Thursday, November 16, 2023 12:49 PM

To: Bard, Joel < Joel.Bard@pfizer.com >

Cc: Jackobek, Ryan Andrew <a href="Ryan.Jackobek@pfizer.com">Ryan.Jackobek@pfizer.com</a>; Francis, Christopher Alan <a href="Christopher.Francis@pfizer.com">Christopher.Francis@pfizer.com</a>; Bitzas, Gabrielle <a href="Gabrielle.Bitzas@pfizer.com">Gabrielle.Bitzas@pfizer.com</a>; Francis, Christopher Alan <a href="Christopher.Francis@pfizer.com">Christopher.Francis@pfizer.com</a>; Bitzas, Gabrielle <a href="Gabrielle.Bitzas@pfizer.com">Gabrielle.Bitzas@pfizer.com</a>; Francis, Christopher.Francis@pfizer.com</a>; Bitzas, Gabrielle <a href="Gabrielle.Bitzas@pfizer.com">Gabrielle.Bitzas@pfizer.com</a>; Francis.

Kerns, Kelvin < Kelvin.Kerns@pfizer.com >

Subject: RE: GeneData Error for new VSIG4 sequences

Hi Joel,

May you help with the complexation issue detailed below? This Issue is happening for GBT-VSIG4-1044, 1047, and 1049. Please let me know if you have any further questions on what is going on.

Thanks, Chris

**Christopher Sausen** 

Senior Scientist

BioMedicine Design

**Biotherapeutic Production** 

O: 1-978-247-1374 M: 1-218-851-7391

Christopher.sausen@pfizer.com



From: Jackobek, Ryan Andrew < Ryan. Jackobek@pfizer.com >

Sent: Thursday, November 16, 2023 11:31 AM

 $\textbf{To:} \ \ \textbf{Sausen, Christopher.Sausen@pfizer.com}{>}; \ \textbf{Francis, Christopher.Alan}{<} \underline{\textbf{Christopher.Francis@pfizer.com}}{>}; \ \textbf{Bitzas, Gabrielle.Bitzas@pfizer.com}{>}; \ \textbf{Francis, Christopher.Alan}{<} \underline{\textbf{Christopher.Francis@pfizer.com}}{>}; \ \textbf{Bitzas, Gabrielle.Bitzas@pfizer.com}{>}; \ \textbf{Francis, Christopher.Alan}{<} \underline{\textbf{Christopher.Francis@pfizer.com}}{>}; \ \textbf{Bitzas, Gabrielle.Bitzas@pfizer.com}{>}; \ \textbf{Francis, Christopher.Bitzas.max}{>}; \ \textbf{Francis, Christopher.Bitz$ 

Kerns, Kelvin < <u>Kelvin.Kerns@pfizer.com</u>>

Subject: RE: GeneData Error for new VSIG4 sequences

Hi Chris,

Complexation needs to be changed to 2, but it looks like its grayed out. Maybe a Joel fix.

Ryan

From: Sausen, Christopher < <a href="mailto:Christopher.Sausen@pfizer.com">Christopher.Sausen@pfizer.com</a>

Sent: Thursday, November 16, 2023 11:28 AM

To: Francis, Christopher Alan < <a href="mailto:Christopher.Francis@pfizer.com">"> Bitzas, Gabrielle < <a href="mailto:Gabrielle.Bitzas@pfizer.com">"> Jackobek, Ryan Andrew < <a href="mailto:Ryan.Jackobek@pfizer.com">Ryan.Jackobek@pfizer.com</a> ; Bitzas, Gabrielle < <a href="mailto:Gabrielle.Bitzas@pfizer.com">Gabrielle.Bitzas@pfizer.com</a> ; Jackobek, Ryan Andrew < <a href="mailto:Ryan.Jackobek@pfizer.com">Ryan.Jackobek@pfizer.com</a> ; Bitzas, Gabrielle < <a href="mailto:Gabrielle.Bitzas@pfizer.com">Gabrielle.Bitzas@pfizer.com</a> ; Jackobek, Ryan Andrew < <a href="mailto:Ryan.Jackobek@pfizer.com">Ryan.Jackobek@pfizer.com</a> ; Bitzas, Gabrielle < <a href="mailto:Gabrielle.Bitzas@pfizer.com">Gabrielle.Bitzas@pfizer.com</a> ; Jackobek, Ryan Andrew < <a href="mailto:Ryan.Jackobek@pfizer.com">Ryan.Jackobek@pfizer.com</a> ; Bitzas, Gabrielle < <a href="mailto:Gabrielle.Bitzas@pfizer.com">Gabrielle.Bitzas@pfizer.com</a> ; Jackobek, Ryan Andrew < <a href="mailto:Ryan.Jackobek@pfizer.com">Ryan.Jackobek@pfizer.com</a> > <a href="mailto:Gabrielle.Bitzas@pfizer.com">Ryan.Jackobek@pfizer.com</a> > <a href="mailto:Gabrielle.

Subject: GeneData Error for new VSIG4 sequences

Hi expression,

For the new VSIG4 designs, in GeneData the complex says it is a 1-mer, and lists data for a 1-mer instead of a 2-mer (see example of 1044 below). Can this be fixed to

show that the complex is a 2-mer?

Thanks, Chris

## Target Product TPP-112645 (GBT-VSIG4-1044) (pfizer.com)



Absorbance [280nm] of 1mg/mL 1.61

Isoelectric Point, pl 6.6

Calculation Method pl from EMBOSS Pepstats version 6.5.7.0a. Mass calculated in GDBxT is oxidized, deglycosylated.

Sequence

1 RPILEVPESV TCPWKGDCNL PCTYDPLQGY TQVLVKWLVQ RGSDPVTIFL RDSSGDHIQQ AKYQGRLHVS HKVPGDVSL(
81 LSTLEMDDRS HYTCEVTWQT PDGNQVVRDK ITELRVQKHS SKLLATATEA PTTMTYPLKA TSTVKQSWDW TTDMDGYLGI
161 TSAGPGASLP GGGGEPKSSD KTHTCPPCPA PELLGGPSVF LFPPKPKDTL MISRTPEVTC VVVDVSHEDP EVKFNWYVD

241 VEVHNAKTKP REEQYNSTYR VVSVLTVLHQ DWLNGKEYKC KVSNKALPAP IEKTISKAKG QPREPQVYTL PPSREEMTKI

321 QVSLTCLVKG FYPSDIAVEW ESNGQPENNY KTTPPVLDSD GSFFLYSKLT VDKSRWQQGN VFSCSVMHEA LHNHYTQKSI

401 SLSPGK

Now empty properties

The complex sequence has been generated in-silico by concatenating the monomer sequence 1 times (complexation number).

# **Christopher Sausen**

Senior Scientist

#### **BioMedicine Design**

**Biotherapeutic Production** 

**O**: 1-978-247-1374 **M**: 1-218-851-7391 <u>Christopher.sausen@pfizer.com</u>

